## Unity® and LAMP Words for Life™ – Choosing a device

Language Acquisition through Motor Planning (LAMP) is a therapeutic approach based on neurological and motor learning principles. LAMP therapy focuses on giving the individual independent access to vocabulary on voice output AAC devices that use consistent motor plans for accessing vocabulary. Vocabulary targets are selected collaboratively with the family, based on the individual's interests and desires. Multisensory input including consistent movement, vision and auditory feedback with natural consequences in their normal environment occurs to maximise language acquisition.

The LAMP approach was first developed by utilising motor learning principles along with the <u>Unity®</u> language system. The <u>LAMP</u> <u>Words for Life™</u> language system was developed to combine the power of Unity with the LAMP principles and strategies. Both language systems strictly utilise the principles of motor planning and can be used in conjunction with the LAMP approach.

Liberator offers several devices with Unity® and/or LAMP Words for Life $^{TM}$ . All Liberator dedicated devices are available for a free 3 week trial with support from our consultants. The LAMP Words for Life $^{TM}$  app is available through the App store.

	Accent 800/1000/1400	Liberator Rugged 8	iPad with LAMP WFL
			Ilike playing  I i i i i i i i i i i i i i i i i i i
Unity	<b>✓</b>		?
LAMP Words for Life	<b>✓</b>	<i>V</i>	<b>✓</b>
Screen size, weight	8/10/14 inch, 1090 - 2680 g	8 inch, 780g	Various
Operating system	Windows (NuVoice)	Android (Chat software)	iOS
84 location keyguard available	V	<i>V</i>	Through external suppliers
Stand, handle & durable case included	V	V	2
Hard-wired enhanced amplification	V	<b>✓</b>	2
Additional access options	Switch Scanning Head tracking Eyegaze (1000 & 1400)	2	2
Choice of dedicated or integrated device	V	<i>V</i>	2
Menu locking	Password or USB override	Password	Guided Access
Data logging	<b>✓</b>	V	?
Direct interface with word processing programs	V	Limited	2
Warranty	3 year	3 year Includes accidental damage	Options through Apple or external agencies
Free Liberator consultant support	V	<b>✓</b>	2
	Backup to internal storage	Backup to internal storage	Backup via iTunes or AirDrop
	Vocabulary transfer to computer via USB or cable	Vocabulary transfer to computer via cable	
Additional features	Editing & resource development using free Windows support software (NuVoice PASS)	Editing & resource development using free Windows support software (Chat Editor)	
	Computer keyboard and mouse emulation	Water resistant, proof rating of IP 66, tempered glass	
	Infrared environmental controls		

