

Nova Chat Vocabulary Options Reference Guide

VOCABULARY	BUTTON BOX LAYOUT	AGE GROUP	TARGET AUDIENCE	VOCABULARY FEATURES
WordPower20 Simply	20	Child Adolescent Adult	-motor problems -vision issues -cognitive delays	<ul style="list-style-type: none"> • Phrase based • Word based • Symbols located on most buttons • Different morphological endings
WordPower24	24	Child Adolescent Adult	-broad spectrum of needs	<ul style="list-style-type: none"> • Logical next words • Very Dynamic • Different morphological endings
WordPower24 with Phrases	24	Child Adolescent Adult	- motor issues -cognitive delays	<ul style="list-style-type: none"> • Want to use when word based isn't working • Phrase based
WordPower30	30	Child Adolescent Adult	-different access methods	<ul style="list-style-type: none"> • Designed for key guard and scanning
WordPower42	42	Child Adolescent Adult	-broad spectrum of needs	<ul style="list-style-type: none"> • Very dynamic • Very predictive
WordPower42 Basic	42	Child Adolescent Adult	-broad spectrum of needs	<ul style="list-style-type: none"> • More simplified version of ChatPower 42 • More static screens • More verbs on the home page (less transitions)
WordPower48	48	Child Adolescent Adult	-different access methods	<ul style="list-style-type: none"> • Designed for key guard and scanning
WordPower80	80	Child Adolescent Adult	-broad spectrum of needs	<ul style="list-style-type: none"> • Good for those who like to spell • More core words

Nova Chat Vocabulary Options Reference Guide

myQuickChat	4 location 8 location 12 location	Child Adult	-complex communication needs	<ul style="list-style-type: none"> • Each vocab has an identical 16 location home page • Phrase based
VocabPC	12 location	Adolescents Adults	-developmental disabilities	<ul style="list-style-type: none"> • Carrier phrases • Interactive sentences • Activity vocabulary • Naming words
MultiChat 15 Student	15 Location	School aged children	-emerging language skills	<ul style="list-style-type: none"> • Sentences • Phrases • Individual words • Visual scenes • News-2 -You starter page
MultiChat 15 Adolescent	15 Location	Adolescents	-emerging language skills	<ul style="list-style-type: none"> • Sentences • Phrases • Individual words • Visual scenes • Pages for work, shopping & socialization
MultiChat 15 Adult	15 Location	Adults	-emerging language skills	<ul style="list-style-type: none"> • Sentences • Phrases • Individual words • Visual scenes • Pages for work and socialization
4-Basic	4 Location		-beginning communicators -used for evaluations	
Spelling	Keyboard page			<ul style="list-style-type: none"> • Four word prediction buttons • Some pre-stored phrases