Communication **Device Guide**



Accent® 800

Accent® 1000

Accent® 1400

Liberator Rugged 10









Android OS

8 inches / 20 cm

Liberator Rugged 8



Android OS

10 inches / 26 cm

Coffware	

NuVoice

Please see the Software Guide for vocabularies available

Operating System	Windows	Windows	Windov	
Screen size	8 inches / 20 cm	10 inches / 26 cm	14 inches /	
Screen durability	Gorilla glass	Gorilla glass	Heat Tempere	
Weight	1.13 kg	1.5 kg	2.7 kg	
USB ports	3 x USB, 1 x USB-C	3 x USB, 1 x USB-C	3 x US	
Keyguards / Touchguides	Keyguards & Touchguides	Keyguards & Touchguides Also 84 Tactile Keyguard	Keyguards & To	
Eye Gaze	No	Yes – Look 1	Yes – Loc	
Head Pointing	Not recommended	Yes – NuPoint	Yes – NuF	
Switch Ports	Yes - 2	`Yes - 2	Yes - 2	
Environmental Control	Infrared + in-built pages that work with voice assistants (e.g. Alexa).			
Mounting	Compatible with Quick Release (Rehadapt) and DAESSY mounting systems			
Battery Life	9-10 hours	9-10 hours 6-8 hours running LOOK	13-15 ho 10-12 hours run	

	Windows			
cm	10 inches / 26 cm			
s	Gorilla glass			
	1.5 kg			
SB-C	3 x USB, 1 x USB-C			
hguides	Keyguards & Touchguides Also 84 Tactile Keyguard			
	Yes – Look 1			
ded	Yes – NuPoint			
	`Yes - 2			
ared + in-built pages that work with voice assis				

	Windows					
	14 inches / 36 cm					
	Heat Tempered Glass					
	2.7 kg					
	3 x USB					
	Keyguards & Touchguides					
	Yes – Look 2					
	Yes – NuPoint					
	Yes - 2					
istar	istants (e.g. Alexa).					

Gorilla glass	Gorilla glass			
840 grams	1.23 kg			
1 x USB-C	1 x USB-C			
Standard keyguards available for certain grid sizes. Custom keyguards also available.				
No	No			
No	No			
No	No			
Pages available that work wi	ith voice assistants (e.g., Alexa).			

Chat Software - LAMP Words for Life vocabularies

and/or Chat vocabularies Please see the Software Guide for vocabularies available

Up to 15 hours with typical AAC use

Up to 15 hours with typical AAC use

Accessories

10-12 hours running LOOK Additional accessories are available. Please visit https://liberator.net.au/store/

13-15 hours

Communication Device Guide



NovaChat® 8

NovaChat® 10

VersaEye Bundle

Via Mini









Software	Chat Software Please see the Software Guide for vocabularies available		iOS Apps Please see the Software Guide for vocabularies available	iOS Apps Please see the Software Guide for vocabularies available
Operating System	Android OS	Android OS	IOS	IOS
Screen size	8 inches / 20 cm	11 inches / 28 cm	12.9 inches / 33 cm	8.3 inch / 21 cm
Screen durability	Gorilla Glass	Gorilla Glass	Aluminosilicate glass	Aluminosilicate glass
Weight	950 grams	1.25 kg	1.5 kg with iPad Pro 12.9"	816 grams
USB ports	Yes – 1 (USB-C)	Yes – 1 (USB-C)	Yes – 1 (USB-C)	Yes – 1 (USB-C)
Keyguards / Touchguides	, 5	lable for certain grid sizes. Duchguides also available.	Yes, through alternate providers.	Standard keyguards available for certain grid sizes.
Eye Gaze	No	No	Hiru eye tracker	No
Head Pointing	No	No	Yes - through iPad Pro	Yes - through iPad Mini
Switch Ports	Yes – 2	Yes – 2	Yes – 1	Yes – 2
Environmental Control	Pages available that work with voice assistants (e.g., Alexa).		Pages available that work with voice assistants (e.g., Alexa).	Pages available that work with voice assistants (e.g. Alexa)
Mounting	Compatible with Quick Release (Rehadapt) and DAESSY mounting systems		Versa QRM Mounting Plate	QRM Mounting Plate
Battery Life	8 hours 45 mins (with continuous use)	10 hours (with continuous use)	9 hours typical battery life	10 hours typical battery life
Accessories	Additional accessories are available. Please visit https://liberator.net.au/store/			

Software Guide (Vocabularies)

NuVoice vocabularies

- Unity: 28, 36, 45, 60, 84, 144 Sequenced or One Hit
- LAMP Word for Life (WFL): 84 Full, Transition and One Hit
- LAMP WFL Vision Impairment (VI): 84 Full, Transition or One Hit
- Essence: 36, 60, 84, 144 ABC or QWERTY
- CoreScanner: 1, 8, 23, 25, 108. 190, 1000s core words
- WordPower for Accent: Word or Pictures 28, 36, 45, 60, 84 (SS or PCS)
- WordCore: 45, 84, 84 Scan, 144)
- Page-Based: 8 Phrase, 8 Auditory Scan, 15 Talk, 45 ABC with Pages

Chat software - LAMP Words for Life

• LAMP WFL 84 Full. Transition and One Hit

Chat software - Chat vocabularies

SymbolStix (SS) or Picture Communication Symbols $^{\circledR}$ (PSC) available

- WordPower: 20, 25, 42 Basic, 42, 48, 60 Basic, 60, 80, 108, 108 w keyboard, 140 Scan & Touch
- VocabPC Adolescent-Adult 12
- Spelling
- myQuickChat: Child or Adult 4, 8, 12
- myCore 117
- Multichat 15 (Student, Adolescent, or Adult)
- Communication Journey: Aphasia 12
- 4 Basic
- Basic Scan 8
- Intermediate Scan 8

iOS Apps

Via Mini and VersaEye bundle comes with one of the following vocabulary apps

- LAMP Words for Life
- TouchChat HD- AAC w/ WordPower
- Dialogue AAC

Additional accessories are available. Please visit https://liberator.net.au/store/