# **Communication Device Guide**



Accent <sup>®</sup> 800	Accent <sup>®</sup> 1000	Accent <sup>®</sup> 1400
Hola	Hele	\$0144 A
		MORE MORE Events Events Events Events MORE
		A S D F G H J K

SOFTWARE	<b>NUVOICE</b> Please see the Software Guide for vocabularies available					
Operating System	Windows	Windows	Windows			
Screen size	8 inches / 20 cm	10 inches / 26 cm	14 inches / 36 cm			
Screen durability	Gorilla glass	Gorilla glass	Heat Tempered Glass			
Weight	1.13 kg	1.5 kg	2.7 kg			
JSB ports	3 x USB, 1 x USB-C	3 x USB, 1 x USB-C	3 x USB			
Keyguards / Touchguides	Keyguards & Touchguides	Keyguards & Touchguides Also 84 Tactile Keyguide	Keyguards & Touchguides			
Eye Gaze	No	Yes – Look 1	Yes – Look 2 Yes – NuPoint			
lead Pointing	Not recommended	Yes – NuPoint				
Switch Ports	Yes - 2	Yes - 2	Yes - 2			
Environmental Control	Infrared + in-built pages that work with voice assistants (e.g. Alexa).					
Mounting	Compatible with Quick Release (Rehadapt) and DAESSY mounting systems					
Battery Life	9-10 hours	9-10 hours 6-8 hours running LOOK	13-15 hours 10-12 hours running LOOK			
Accessories	Additional accessories are available. Please visit https://liberator.net.au/store/					

#### SOFTWARE GUIDE (Vocabularies)

NuVoice vocabularies

- Unity: 28, 36, 45, 60, 84, 144 Sequenced or 1-Hit
- LAMP Word for Life™ (WFL): 84 Full, Transition and 1-Hit
- LAMP WFL Vision Impairment (VI): 84 Full, Transition or 1-Hit
- Essence: 36, 60, 84, 144 ABC or QWERTY
- CoreScanner: 1, 8, 23, 25, 108. 190, 1000s core words
- WordPower for Accent: Word or Pictures 28, 36, 45, 60, 84 (SS or PCS)
- WordCore: 45, 84, 84 Scan, 144)
- Page-Based: 8 Phrase, 8 Auditory Scan, 15 Talk, 45 ABC with Pages



Not sure which device to try? Talk to one of our Consultants: liberator.net.au/contact-us | (02) 9124 9945

	NovaRugged <sup>®</sup> 8	NovaRugged <sup>®</sup> 10	VersaEye Bundle	Via Mini	
SOFTWARE	<b>CHAT SOFTWARE</b> Please see the Software Guide for vocabularies available	<b>CHAT SOFTWARE</b> Please see the Software Guide for vocabularies available	<b>IOS APPS</b> Please see the Software Guide for vocabularies available	<b>IOS APPS</b> Please see the Software Guide for vocabularies available	
Operating System	Android OS	Android OS	IOS	IOS	
Screen size	8 inches / 20 cm	11 inches / 28 cm	12.9 inches / 33 cm	8.3 inch / 21 cm	
Screen durability	Gorilla Glass	Gorilla Glass	Aluminosilicate glass	Aluminosilicate glass	
Weight	952 grams	1.25 kg	1.5 kg with iPad Pro 12.9"	816 grams	
USB ports	1 x USB-C	1 x USB-C	Yes – 1 (USB-C)	Yes – 1 (USB-C)	
Keyguards / Touchguides	Standard keyguards available for certain grid sizes	Standard keyguards available for certain grid sizes	Yes, through alternate providers.	Standard keyguards available for certain grid sizes.	
Eye Gaze	No	No	Hiru eye tracker	No	
Head Pointing	No	No	Yes - through iPad Pro	Yes - through iPad Mini	
Switch Ports	2 Points	2 Points	Yes – 1	Yes – 2	
Environmental Control	No		Pages available that work with voice assistants (e.g., Alexa).	Pages available that work with voice assistants (e.g. Alexa)	
Mounting	Compatible with Quick Release (Rehadapt) and DAESSY mounting systems		Versa QRM Mounting Plate	QRM Mounting Plate	
Battery Life	Up to 12 hours with typical AAC use, 8 hours continuous	Up to 12 hours with typical AAC use, 10 hours continuous	9 hours typical battery life	10 hours typical battery life	
Accessories	Additional accessories are available. Please visit https://liberator.net.au/store/				

### **CHAT SOFTWARE**

#### Lamp Words for Life<sup>™</sup> vocabularies

- LAMP WFL Full
- LAMP WFL Transition
- LAMP WFL 1 Hit

#### **Chat Vocabularies**

- SymbolStix (SS) or Picture Communication Symbols® (PCS) available
- WordPower: 20, 25, 42 Basic, 42, 48, 60 Basic, 60, 80, 108, 108 w keyboard, 140 Scan & Touch
- VocabPC Adolescent-Adult 12
- Spelling
- myQuickChat: Child or Adult 4, 8, 12
- myCore 117 -
- Multichat 15 (Student, Adolescent, or Adult) -
- Communication Journey: Aphasia 12
- 4 Basic -
- Basic Scan 8
- Intermediate Scan 8

#### **iOS** Apps

Via Mini and VersaEye bundle comes with one of the following vocabulary apps

- LAMP Words for Life<sup>™</sup>
- TouchChat HD- AAC w/ WordPower
- Dialogue AAC



## Got questions? Book a chat with one of our Consultants: liberator.net.au/contact-us | (02) 9124 9945

© 2025 PRC-Saltillo. Non-commercial reprint rights for clinical or personal use granted with inclusion of copyright notice. Commercial use prohibited; may not be used for resale. WordPower® is a trademark of Inman Innovations, Inc ©2001-2025. PCS® is a trademark of Tobii Dynavox LLC. The symbols used in the work are owned by n2y, LLC. Used with permission under license. ©2024. Product prices and features are subject to change without notice.